

**Concentration PPM (mg/Liter) vs %Transmittance  
0-1,500 PPM COD Vials**

**620 nm**

%T	%T UNITS									
TENS	0	1	2	3	4	5	6	7	8	9
10										1657
20	1605	1557	1510	1466	1423	1382	1343	1305	1269	1234
30	1200	1167	1136	1105	1075	1046	1018	991	964	938
40	913	888	864	840	817	795	773	752	730	710
50	690	670	650	631	613	594	576	559	541	524
60	507	491	475	459	443	427	412	397	382	368
70	353	339	325	311	298	284	271	258	245	233
80	220	208	195	183	171	159	148	136	125	113
90	102	91	80	70	59	48	38	27	17	7

**PPM (mg/Liter) = 2301(abs) -3**

**CHEMetrics, Incorporated  
Route 28, Calverton, VA 20138  
Phone (540) 788-9026  
www.chemetrics.com**

Oct. 04, Rev. 5