

**Concentration vs %Transmittance**  
**Chlorine 1, PPM (mg/Liter)**

Cat. No. K-2503

565 nm

TENS	%T									
	0	1	2	3	4	5	6	7	8	9
10	3.49	3.35	3.22	3.09	2.98	2.88	2.78	2.69	2.60	2.52
20	2.45	2.37	2.30	2.23	2.17	2.11	2.05	1.99	1.94	1.89
30	1.83	1.78	1.74	1.69	1.65	1.60	1.56	1.52	1.48	1.44
40	1.40	1.36	1.33	1.29	1.26	1.22	1.19	1.16	1.13	1.10
50	1.06	1.03	1.01	0.98	0.95	0.92	0.89	0.87	0.84	0.82
60	0.79	0.76	0.74	0.72	0.69	0.67	0.65	0.62	0.60	0.58
70	0.56	0.54	0.52	0.49	0.47	0.45	0.43	0.41	0.39	0.38
80	0.36									
90										

PPM (mg/Liter) = 3.47(abs) +0.02

CHEMetrics, Incorporated  
 Route 28, Calverton, VA 20138  
 Phone (540) 788-9026  
 www.chemetrics.com

June 03, Rev. 10